


Graphing



Item No.	Qty	Description	Inventory Code	Material
D1	1	1L SU Storage Bags	ACAX0012	ULDPE
B1	2	1/4" 3/8" (6.4*9.6mm) L-Flex Tubing, Length 500mm	DBBX0063	TPE
B2	1	1/4" 3/8" (6.4*9.6mm) Silicone Tubing (LePure), Length 200mm	DAAX0066	Silicone
C1	1	1/4" (6.4mm) Hose Barb Non-Valved in-line coupling body	EBAX0005	PC
C2	1	MPC coupling plug	EBAX0028	PC
C3	1	1/4" (6.4mm) Hose barb non-valved in-line coupling insert	EBAX0017	PC
C4	1	MPC coupling cap	EBAX0025	PC
C5	1	1/4" (6.4mm) Female luer thread barb	EAAAX0143	PP
C6	1	needle-less ampler	EABX0040	PC
C7	1	needle-less sampler cap	EABX0005	PC
F1	6	2.5mm cable tie	ADAA0005	Nylon
F2	3	3.6mm cable tie	ADAA0004	Nylon
F3	3	Robert clip (small)	EDBA0026	PA

Customer Approval		LePure Biotech Co., Ltd.				
Company:		Drawing No. HP001-4AK				
Signature:		Title: SU Storage Bags				
Date:		Designed By:	Date:	Checked By:	Date:	Edition
		Jiajia Guo	2023.03.12	Alan	2023.03.12	Rev01
<small>Design Approval the manufacturing of the system(s) based on this design. By signing, you also assume full</small>		<small>We request that you do not disclose this information to any third party involved in this drawing, in order to protect our customer's confidentiality and intellectual property</small>				