

## ABCB7 Polyclonal Antibody

Catalog number: 11158-1-AP


Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

ABCB7; ABC transporter 7 protein, ABC7, ABCB7, ASAT, Atm1p, EST140535



mouse liver tissue were subjected to SDS PAGE followed by western blot with 11158-1-AP(ABCB7 antibody) at dilution of 1:1000



IP Result of anti-ABCB7 (IP:11158-1-AP, 4µg; Detection:11158-1-AP 1:500) with mouse liver tissue lysate 4000µg.

### Applications

Tested applications:	ELISA, WB, IP
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated ABCB7 MW:	753aa,83 kDa
Observed ABCB7 MW:	83kd
Positive WB detected in	Mouse liver tissue, mouse thymus tissue
Positive IP detected in	Mouse liver tissue
Recommended dilution:	WB: 1:200-1:2000 IP: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag1544
GenBank accession number:	BC006323
Gene ID (NCBI):	22
Full name:	ATP-binding cassette, sub-family B (MDR/TAP), member 7

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.