

AF9 Polyclonal Antibody

Catalog number: 12825-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

MLLT3; AF9, FLJ2035, MLLT3,

Protein AF 9, YEATS3



Jurkat cells were subjected to SDS PAGE followed by western blot with 12825-1-AP (MLLT3 antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated AF9 MW:	568aa, 63 kDa
Observed AF9 MW:	63kd
Positive WB detected in	Jurkat cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3528
GenBank accession number:	BC036089
Gene ID (NCBI):	4300
Full name:	Myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.