

YOUR DATA, ANYWHERE, ANYTIME




LogTag[®]

LogTag Online
Fixed-Site Cloud Storage Service

LogTag Online (LTO) is a Cloud-based online environmental management system providing alarm notifications, compliance, record storage and reporting functions.

LTO allows its users to download data from LogTag dataloggers, view the data in the cloud, and allows the sharing and securing of the data using a permissions system.



WiFi Compatibility:


Works with all LogTag loggers



UTRED30-WiFi
Not Included



LTI-WiFi
Interface cradle
Not Included



LTI-WM-WiFi
Interface cradle
Not Included

Features



Immediate notification of a temperature excursion sent to your email or phone.



Data is uploaded to the cloud in real time.



Suitable for all users from small clinics to the largest hospitals.



LogTag Online automates vaccine temperature monitoring and maintains compliance to VFC and CDC Guidelines.




Access your data through LogTag Online on any device with an internet connection and internet browser.



Setting up your account is quick and easy, then you can access your data online.



LogTag Online can be accessed using any LogTag data logger, so you don't need to upgrade your existing loggers in order to get started.



Support for LogTag Online is available through LogTag's Global distribution network.


Applications



Vaccines / Medical & Healthcare



Industrial & Manufacturing



Warehouse



Food & Produce



Storage Rooms



Logistics

Account Structure

Feature	Basic (Free) Account (No Activation Code)	Standard (Paid) Account* (Activation Code required)
Number of Locations	1 Limited Location	1 Full Location per purchased Activation Code
Location History	3 Months	10 Years
Notifications	Email Only	Email + SMS
SMS (USA)	None	Unlimited
SMS (Outside USA)	None	Up to 25 messages per month
Maximum number of files	30	1000
File Storage	1 Year	10 Years
Number of registered devices	Unlimited	Unlimited
Number of users	2	Unlimited

*To activate and use LogTag Online features (depending on user requirement), contact your local provider to purchase activation code(s) – and be on a Standard (Paid) account.

Security Compliance

LogTag Online is committed to providing you with a trusted set of cloud services. LogTag Online complies with all requirements for electronic records including [FDA 21 CFR Part 11](#). We have leveraged our decades-long industry experience building enterprise to create a robust set of LogTag Online security technologies and practices. These work to help reduce the cost, complexity, and risk associated with security in the cloud.

To download the full LogTag Online Security Compliance PDF, go to the [LTO Compliance tab](#).