

## HTRA4 Polyclonal Antibody

Catalog number: 14095-1-AP


Size: 33 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

HTRA4; FLJ90724



human placenta tissue were subjected to SDS PAGE followed by western blot with 14095-1-AP(HTRA4 Antibody) at dilution of 1:600

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human placenta tissue
Calculated HTRA4 MW:	476aa,51kd
Observed HTRA4 MW:	51kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag5227
GenBank accession number:	BC057765
Gene ID (NCBI):	203100
Full name:	HtrA serine peptidase 4

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.