

ITPK1 Polyclonal Antibody

Catalog number: 10568-1-AP

Size: 39 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

ITPK1; Ins(1,3,4)P(3) 5/6 kinase,

ITPK1, ITRPK1



human brain tissue were subjected to SDS PAGE followed by western blot with 10568-1-AP(ITPK1 antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated ITPK1 MW:	48 kDa
Observed ITPK1 MW:	45 kDa
Positive WB detected in	Human brain tissue, HEK-293 cells, mouse eye tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag0834
GenBank accession number:	BC007428
Gene ID (NCBI):	3705
Full name:	Inositol 1,3,4-triphosphate 5/6 kinase

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.