



The Adtemp™ 412

Thomas No. 20A00H120



Our most cost effective digital stick thermometer, now with beeper.

- Features:**
- Super low-cost digital thermometer
 - Compact size: less than 5" long, weighs just 14 grams
 - Range 90.0°F to 109.9°F \pm .2°F or 32.0°C to 43.9°C \pm .1°C depending upon scale selection
 - Non-replaceable 1.55V (LR41) type battery provides up to 4,500 measurements
 - Integral storage case keeps instrument sanitary
 - Bulk unit packaged 20/box, 240/case
 - Meets the accuracy requirements for intermittent thermometers under ASTM E1112

The Adtemp™ Ultra 417

Thomas No. 20A00H121



The ultra-fast stick-type digital thermometer uses predictive technology for a two-second oral reading.

- Features:**
- Two-second measurement using proprietary predictive technology
 - Range 95°F to 107.6°F \pm .2°F or 35°C to 42°C \pm .1°C depending upon scale selection
 - Replaceable 3V (2032Li) coin-type battery provides a minimum of 1,000 full-duty measurements
 - Integral scabbard-type carry case
 - Includes five probe sheaths
 - Contemporary Euro design

The Adtemp™ 415 Flex

Thomas No. 20A00H122



Quick read digital stick thermometer w/ flexible tip for patient safety.

- Features:**
- 18 second oral, 10 second rectal temperature measurement
 - Large LCD for easy reading
 - Range 90°F to 109.9°F \pm .2°F or 32°C to 43.9°C \pm .1°C depending upon scale selection
 - Flexible probe tip
 - Replaceable 1.55v (LR41) type battery provides up to 4,000 measurements
 - Includes 5 probe sheaths
 - Clear plastic scabbard storage case
 - Each in box, 12/carton, 144/case
 - Meets the accuracy requirements for intermittent thermometers under ASTM E1112

Note: These products are not available for purchase in Canada.

Disclaimer: Pricing and offers subject to change without notice.