

RARA Polyclonal Antibody

Catalog number: 10331-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

RARA; NR1B1; RAR

Applications

Tested applications: ELISA
Species specificity: Human, mouse, rat; other species not tested.
Calculated RARA MW: 51kd
Observed RARA MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen: Ag0420
GenBank accession number: BC008727
Gene ID (NCBI): 5914
Full name: Retinoic acid receptor, alpha

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.