

TEKT2 Polyclonal Antibody

Catalog number: 13518-1-AP


Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

TEKT2; TEKTB1; TEKTIN-T; h-
tektin-t



mouse brain tissue were subjected to SDS PAGE followed by western blot with 13518-1-AP (TEKT2 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human trachea using 13518-1-AP (TEKT2 Antibody) at Dilution 1:50 (under 40x lens)

Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse brain tissue
Calculated TEKT2 MW:	430aa, 50kd
Observed TEKT2 MW:	54kd
Positive IHC detected in:	Human trachea
Recommended dilution:	WB: 1:200-1:2000 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4457
GenBank accession number:	BC035620
Gene ID (NCBI):	27285
Full name:	Tektin 2 (testicular)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.