

Levo Pipette Controller





-----Compress the suction bellows before attaching the pipette

-----Button to expel residual liquid when using blow-out pipets

-----Filling and delivery are easily controlled with the small sensitive lever

Ordering information

Thomas No.	Cat No.	Description
1227F67	728001	Levo - Blue
1227F68	728002	Levo - Red
1227F69	728003	Levo - Green
1227F70	728004	Levo - Yellow
1227F71	728005	Levo - Pink

Features:

- Comfortable and simple to use
- Precise pipetting control
- Robust and lightweight
- 1-100ml volumetric or serological pipettes can be used
- Uses integral and replaceable 3um filter
- Easy to maintain and clean
- Five colors
- Autoclavable silicone pipette holder

DRAGONLAB