

## SSX2IP Polyclonal Antibody

Catalog number: 13694-1-AP

Size: 44 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

SSX2IP; ADIP; FLJ10848;

KIAA0923; MGC75026



HEK-293 cells were subjected to SDS PAGE followed by western blot with 13694-1-AP(SSX2IP Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	HEK-293 cells
Calculated SSX2IP MW:	614aa,71kd
Observed SSX2IP MW:	64kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag4584
GenBank accession number:	BC033637
Gene ID (NCBI):	117178
Full name:	Synovial sarcoma, X breakpoint 2 interacting protein

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.