

NXT2 Polyclonal Antibody

Catalog number: 12136-1-AP


Size: 24 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

NXT2; DC9, NTF2 related
export protein 2, NXT2, P15 2,
Protein p15 2



Immunohistochemical of paraffin-
embedded human colon cancer using
12136-1-AP(NXT2 antibody) at
dilution of 1:50 (under 10x lens)

Background

NXT2 is a member of NXT family proteins that are generally involved in exporting nuclear RNA in eukaryotic cells. The function of NXT2 has not been widely studied, and is yet to be fully elucidated.

Applications

Tested applications:	ELISA, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated NXT2 MW:	197aa, 23 kDa
Observed NXT2 MW:	
Positive IHC detected in	Human colon cancer tissue
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2782
GenBank accession number:	BC014888
Gene ID (NCBI):	55916
Full name:	Nuclear transport factor 2-like export factor 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.