

USP25 Polyclonal Antibody

Catalog number: 12199-1-AP

Size: 26 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

USP25; Deubiquitinating

enzyme 25, Ubiquitin

thiolesterase 25, USP on

chromosome 21, USP21, USP25



BxPC-3 cells were subjected to SDS PAGE followed by western blot with 12199-1-AP(USP25 antibody) at dilution of 1:400

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated USP25 MW:	125 kDa
Observed USP25 MW:	125 kDa
Positive WB detected in	BxPC-3 cells, human cerebellum tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2837
GenBank accession number:	BC015930
Gene ID (NCBI):	29761
Full name:	Ubiquitin specific peptidase 25

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.