

GRK3 Polyclonal Antibody

Catalog number: 13727-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

ADRBK2; ADRBK2, BARK2,

Beta ARK 2, GRK3



K-562 cells were subjected to SDS PAGE followed by western blot with 13727-1-AP (ADRBK2 antibody) at dilution of 1:400

Applications

| | |
|--------------------------|--|
| Tested applications: | ELISA, WB |
| Species specificity: | Human, Mouse, Rat; other species not tested. |
| Calculated GRK3 MW: | 387aa, 45kd |
| Observed GRK3 MW: | 70kd and 90kd |
| Positive WB detected in: | K-562 cells |
| Recommended dilution: | WB: 1:400 - 1:1600 |

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

| | |
|---------------------------|-------------------------------------|
| Immunogen: | Ag4685 |
| GenBank accession number: | BC036797 |
| Gene ID (NCBI): | 157 |
| Full name: | Adrenergic, beta, receptor kinase 2 |

Product information

| | |
|----------------------|--|
| Purification method: | Antigen affinity purification |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C. |