



Item No.	Qty	Description	Inventory Code	Material
D1	1	500mL SU Storage Bags	ACAX0012	ULDPE
B1	2	1/4"3/8" (6.4*9.6mm) L-Flex Tubing, Length 300mm	DBBX0003	TPE
C1	1	1/4" (6.4mm) Hose barb non-valved in-line coupling body	EAXX0005	PC
C2	1	MPC coupling plug	EAXX0028	PC
C3	1	1/4" (6.4mm) Hose barb non-valved in-line coupling insert	EAXX0017	PC
C4	1	MPC coupling cap	EAXX0025	PC
F1	4	2.5mm cable tie	ADAA0005	Nylon
F2	2	3.6mm cable tie	ADAA0004	Nylon
F3	2	Robert clip (small)	EDBA0026	PA

Customer Approval		LePure Biotech Co., Ltd.				
Company:		Drawing No.	HP0500-3AK			
Signature:		Title:	SU Storage Bags			
Date:		Designed By:	Date:	Checked By:	Date:	
		Jiajia Guo	2023.03.12	Alan	2023.03.12	
					Edition	
					Rev01	