

XPNPEP1 Polyclonal Antibody

Catalog number: 10661-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

XPNPEP1; SAMP; XPNPEP;

XPNPEPL; XPNPEPL1

Applications

Tested applications: ELISA

Species specificity: Human, mouse, rat; other species not tested.

Calculated XPNPEP1 MW: 70kd

Observed XPNPEP1 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen: Ag1015

GenBank accession number: BC007579

Gene ID (NCBI): 7511

Full name: X-prolyl aminopeptidase (aminopeptidase P) 1, soluble

Product information

Purification method: Antigen affinity purification

Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.