

## HMGN4 Polyclonal Antibody

Catalog number: 11686-1-AP

Size: 56 µg/150 µl


Source: Rabbit

Isotype: IgG


Synonyms:

HMGN4; HMG17L3, HMGN4,

NHC



HeLa cells were subjected to SDS PAGE followed by western blot with 11686-1-AP(HMGN4 antibody) at dilution of 1:400



Immunofluorescent analysis of ( 10% Formaldehyde ) fixed HeLa cells using 11686-1-AP( HMGN4 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

### Applications

Tested applications:	ELISA, WB, IF
Species specificity:	Human; other species not tested.
Calculated HMGN4 MW:	90aa,10 kDa
Observed HMGN4 MW:	15 kDa
Positive WB detected in	HeLa cells
Positive IF detected in	HeLa cells
Recommended dilution:	WB: 1:500-1:5000 IF: 1:50-1:500

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag2298
GenBank accession number:	BC001282
Gene ID (NCBI):	10473
Full name:	High mobility group nucleosomal binding domain 4

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.