

## EIF3B Polyclonal Antibody

Catalog number: 10319-1-AP

Size: 82 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

EIF3B; eIF 3 eta, EIF3 ETA, EIF3

P110, EIF3 P116, EIF3B, EIF3S9,

hPrt1, PRT1, Prt1 homolog

### Applications

**Tested applications:** ELISA  
**Species specificity:** Human, Mouse, Rat; other species not tested.  
**Calculated EIF3B MW:** 93 kDa  
**Observed EIF3B MW:**

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

**Immunogen:** Ag0386  
**GenBank accession number:** BC001173  
**Gene ID (NCBI):** 8662  
**Full name:** Eukaryotic translation initiation factor 3, subunit B

### Product information

**Purification method:** Antigen affinity purification  
**Storage:** PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.