




# Puritan medical swabs – meeting healthcare needs for nearly a century.



## WHY PURITAN MEDICAL SWABS?

Puritan is a name known all over the world. We have built our **long-standing reputation** for **innovation** and **quality** by listening to our customers and meeting the evolving needs of the medical fields.













As always, Puritan responds to market needs with high-quality products. We have evolved alongside the healthcare field. Early in our history, the Puritan line was limited, primarily to tongue depressors and cotton tipped applicators. As medicine advanced, we kept pace, expanding our lines and even developing brand-new product categories and swabs for specialized medical applications.

The all-time Puritan best sellers shown above are time-tested examples. Look inside to see our standard product variety. And please remember, your feedback and custom orders are always welcome. The demands of the medical field are and will always be a driving force behind our R&D and our 100% commitment to consistent quality.



## COTTON TIPPED APPLICATORS

We are very particular about our raw materials. For these applicators, we only use US-grown cotton. We coil and card our USP-grade cotton fiber in our factory and inspect every swab to ensure optimum quality and consistency.

|  | Item No.    | Item Description                              | Packaging pack, weight | Handle length x diameter             | Tip length x diameter               | Full Length         |
|--|-------------|---|------------------------|--------------------------------------|-------------------------------------|---------------------|
|    | 803-PC      | Polystyrene handle, standard tip.             | 10/10/100, 9 lbs       | 74.599 x 2.54 mm<br>2.937 x 0.1"     | 17.449 x 5.156 mm<br>0.687 x 0.203" | 76.2 mm<br>3"       |
|  | 25-803 2PC  |   | 10/100/2, 5 lbs        |                                      |                                     |                     |
|    | 803-PC DBL  | Polystyrene handle, double tipped.            | 10/10/50, 7 lbs        | 73.025 x 2.54 mm<br>2.875 x 0.1"     | 17.45 x 4.369 mm<br>0.687 x 0.172"  | 79.375 mm<br>3.125" |
|    | 803-PCL     | Polystyrene handle, large tip.                | 10/5/50, 5 lbs         | 74.599 x 2.54 mm<br>2.937 x 0.1"     | 17.449 x 6.35 mm<br>0.687 x 0.25"   | 76.2 mm<br>3"       |
|    | 803-WC      | Wood handle, standard tip.                    | 10/10/100, 10 lbs      | 73.812 x 2.54 mm<br>2.906 x 0.1"     | 15.875 x 4.775 mm<br>0.625 x 0.188" | 76.2 mm<br>3"       |
|  | 25-803 2WC  |   | 10/100/2, 5 lbs        |                                      |                                     |                     |
|    | 803-WCL     | Wood handle, large tip.                       | 10/5/100, 7 lbs        | 73.812 x 2.54 mm<br>2.906 x 0.1"     | 17.449 x 5.943 mm<br>0.687 x 0.234" | 76.2 mm<br>3"       |
|    | 806-PC      | Polystyrene handle, standard tip.             | 10/10/100, 17 lbs      | 147.624 x 2.54 mm<br>5.812 x 0.1"    | 17.449 x 5.156 mm<br>0.687 x 0.203" | 152.4 mm<br>6"      |
|  | 25-806 1PC  |   | 10/100/1, 6 lbs        |                                      |                                     |                     |
|  | 25-806 2PC  |   | 10/100/2, 7 lbs        |                                      |                                     |                     |
|  | 806-PCL     | Polystyrene handle, large tip.                | 10/5/100, 10 lbs       | 149.225 x 2.54 mm<br>5.875 x 0.1"    | 17.449 x 6.35 mm<br>0.687 x 0.25"   | 152.4 mm<br>6"      |
|  | 806-WC      | Wood handle, standard tip.                    | 10/10/100, 16 lbs      | 147.624 x 2.54 mm<br>5.812 x 0.1"    | 15.875 x 4.775 mm<br>0.625 x 0.188" | 152.4 mm<br>6"      |
|  | 25-806 1WC  |   | 10/100/1, 6 lbs        |                                      |                                     |                     |
|  | 25-806 2WC  |   | 10/100/2, 7 lbs        |                                      |                                     |                     |
|  | 25-806 5WC  |   | 10/20/5, 5 lbs         |                                      |                                     |                     |
|  | 25-806 10WC |   | 10/20/10, 6 lbs        |                                      |                                     |                     |
|  | 806-WCL     | Wood handle, large tip.                       | 10/5/100, 10 lbs       | 147.624 x 2.54 mm<br>5.812 x 0.1"    | 17.449 x 5.943 mm<br>0.687 x 0.234" | 152.4 mm<br>6"      |
|  | 821-WC DBL  | Wood handle, double tipped, tapered one end.  | 10/5/100, 9 lbs        | 147.624 x 2.108 mm<br>5.812 x 0.083" | 15.875 x 4.775 mm<br>0.625 x 0.188" | 152.4 mm<br>6"      |
|  | 823-WC      | Wood handle, tapered mini tip.                | 10/10/100, 8 lbs       | 73.812 x 2.108 mm<br>2.906 x 0.083"  | 12.7 x 3.048 mm<br>0.5 x 0.12"      | 76.2 mm<br>3"       |
|  | 826-WC      | Wood handle, tapered mini tip.                | 10/10/100, 12 lbs      | 147.624 x 2.108 mm<br>5.812 x 0.083" | 12.7 x 3.048 mm<br>0.5 x 0.12"      | 152.4 mm<br>6"      |
|  | 826-WC DBL  | Wood handle, double tipped, tapered mini tip. | 10/10/100, 13 lbs      |                                      |                                     |                     |

continued

## PACKAGING QUANTITY CODES




Our packaging codes work from the outside in. The first number represents the outermost package – this may be a case, a box, or other container. Next is what you'll find when you open the outer packaging – bags, wrappers, tubes, etc. Finally we get to the heart of the packaging – how many swabs are grouped together. Some examples follow:

10/50 = 10 bags of 50 swabs per bag, per case.

10/10/100 = 10 boxes of 10 bags per case, 100 swabs per bag.





10/100/1 = 10 boxes of 100 swabs with each swab individually wrapped or in a tube.

## COTTON TIPPED APPLICATORS, CONTINUED

|   | Item No.   | Item Description                          | Packaging pack, weight | Handle length x diameter             | Tip length x diameter               | Full Length    |
|---|------------|---|------------------------|--------------------------------------|-------------------------------------|----------------|
|  | 830-PC DBL | Paper handle, double tipped.              | 10/10/50, 8 lbs        | 73.025 x 2.54 mm<br>2.875 x 0.1"     | 17.449 x 4.368 mm<br>0.687 x 0.172" | 76.2 mm<br>3"  |
|  | 833-WCS    | Wood handle, small tip, in plastic vials. | 50/100, 8 lbs          | 73.812 x 2.108 mm<br>2.906 x 0.083"  | 15.875 x 3.962 mm<br>0.625 x 0.156" | 76.2 mm<br>3"  |
|  | 869-WCS    | Wood handle, extra small tip.             | 10/10/100, 13 lbs      | 147.624 x 2.108 mm<br>5.812 x 0.083" | 12.7 x 3.175 mm<br>0.5 x 0.125"     | 152.4 mm<br>6" |

## OB/GYN AND PROCTOSCOPIC APPLICATORS








The best swabs – indeed, the best products in any category – are those that are designed to meet the real-world demands of the end user. These swabs are a good example.

|   | Item No.   | Item Description  | Packaging pack, weight | Handle length x diameter            | Tip length x diameter            | Full Length          |
|---|------------|---|------------------------|-------------------------------------|----------------------------------|----------------------|
|  | 808        | OB/GYN applicator, rayon tip, paper handle.               | 10/50, 5 lbs           | 203.2 x 3.962 mm<br>8 x 0.156"      | 28.575 x 12.7 mm<br>1.125 x 0.5" | 207.975 mm<br>8.188" |
|   | 808 100    |   | 10/100, 8 lbs          |                                     |                                  |                      |
|   | 808 BULK   |   | 500 bulk, 4 lbs        |                                     |                                  |                      |
|  | 808 COTTON | OB/GYN applicator, cotton tip, paper handle.              | 10/50, 5 lbs           | 203.2 x 3.962 mm<br>8 x 0.156"      | 28.575 x 12.7 mm<br>1.125 x 0.5" | 207.975 mm<br>8.188" |
|  | 25-808 1PR | OB/GYN applicator, rayon tip, polystyrene handle.         | 10/50/1, 8 lbs         | 196.85 x 3.81 mm<br>7.75 x 0.15"    | 28.575 x 12.7 mm<br>1.125 x 0.5" | 203.2 mm<br>8"       |
|   | 25-808 2PR |   | 10/50/2, 11 lbs        |                                     |                                  |                      |
|  | 816        | Proctoscopic applicator, rayon tip, polypropylene handle. | 10/50, 8 lbs           | 403.225 x 3.81 mm<br>15.875 x 0.15" | 28.575 x 12.7 mm<br>1.125 x 0.5" | 406.4 mm<br>16"      |
|   | 816 100    |   | 10/100, 13 lbs         |                                     |                                  |                      |
|   | 816 BULK   |   | 500 bulk, 6 lbs        |                                     |                                  |                      |

## CYTOLOGY COLLECTION DEVICES

Puritan cytology products are made of the very best materials, in the most effective shapes, to maximize sample collection and provide a comfortable patient experience at the same time.





Products below shown at 80% of actual size.

|  | Item No.   | Item Description   | Packaging pack, weight | Handle length x diameter             | Tip length x diameter               | Full Length          |
|--|------------|--|------------------------|--------------------------------------|-------------------------------------|----------------------|
|    | 2188       | Histobrush, <sup>®</sup> soft nylon bristles, tapered tip. | 10/100, 5 lbs          | 183.362 x 2.971 mm<br>7.219 x 0.117" | 19.05 x 5.08 mm<br>0.75 x 0.2"      | 203.2 mm<br>8"       |
|  | 25-2188    |  | 10/100/1, 8 lbs        |                                      |                                     |                      |
|    | 2194       | Rover's <sup>®</sup> spatula.                              | 15/50, 4 lbs           | 167.5 x 4.4 mm<br>6.594 x 0.173"     | 30 x 16.7 mm<br>1.181 x 0.657"      | 196.85 mm<br>7.75"   |
|    | 2195       | Rover's <sup>®</sup> Cervex-Brush. <sup>®</sup>            | 35/25, 6 lbs           | 166.674 x 4.394 mm<br>6.562 x 0.173" | 28.575 x 19.05 mm<br>1.25 x 0.75"   | 203.2 mm<br>8"       |
|    | 2195 COMBI | Rover's <sup>®</sup> Cervex-Brush <sup>®</sup> combi.      | 4/25, 5 lbs            | 167.335 x 4.394 mm<br>6.588 x 0.173" | 34.29 x 19.05 mm<br>1.35 x 0.75"    | 203.2 mm<br>8"       |
|    | 2196       | Rover's <sup>®</sup> EndoCervex-Brush. <sup>®</sup>        | 10/50, 3 lbs           | 167.487 x 4.394 mm<br>6.594 x 0.173" | 33.782 x 5.435 mm<br>1.333 x 0.214" | 203.2 mm<br>8"       |
|    | 2197       | CerviSoft, <sup>®</sup> soft foam tip.                     | 10/100, 3 lbs          | 177.8 x 2.54 mm<br>7 x 0.1"          | 25.4 x 4.775 mm<br>1 x 0.188"       | 179.374 mm<br>7.062" |
|  | 25-2197    |  | 10/100/1, 7 lbs        |                                      |                                     |                      |
|  | 2199       | Histobrush, <sup>®</sup> soft nylon bristles, tapered tip. | 10/100, 4 lbs          | 157.962 x 3.073 mm<br>6.219 x 0.121" | 19.05 x 5.08 mm<br>0.75 x 0.2"      | 177.8 mm<br>7"       |
|  | 25-2199    |  | 10/100/1, 8 lbs        |                                      |                                     |                      |

## CERVICAL SCRAPERS


Healthcare practitioners have counted on us for years to provide standard wooden scrapers. The quality of the wood, our manufacturing, and quality processes ensure a smooth finish for effective collection and patient comfort.

Products below shown at 30% of actual size.

|  | Item No. | Item Description         | Packaging pack, weight | Full Length    |
|--|----------|--------------------------|------------------------|----------------|
|  | 2173     | Ayres style.             | 10/500, 23 lbs         | 177.8 mm<br>7" |
|  | 25-2173  |                          | 10/100/1, 9 lbs        |                |
|  | 2181     | Medical packaging style. | 10/500, 19 lbs         | 177.8 mm<br>7" |
|  | 2183     | Double-ended.            | 10/500, 21 lbs         | 177.8 mm<br>7" |
|  | 2184     | Lerner style.            | 10/500, 24 lbs         | 177.8 mm<br>7" |

## FOAM WOUND MEASUREMENT DEVICES

The soft foam-tipped applicator is flexible to conform to wound shape and features a metric rule for easy measuring. Perfect for forensics, emergency rooms, long-term care, and more. Choose sterile or non-sterile.

|   | Item No.       | Item Description   | Packaging pack, weight | Handle length x diameter          | Tip length x diameter          | Full Length    |
|---|----------------|--|------------------------|-----------------------------------|--------------------------------|----------------|
|  | 1506-PFB DM    | DM stick, polystyrene handle with printed graduations, foam tip. | 20/50, 2 lbs           | 149.225 x 2.54 mm<br>5.875 x 0.1" | 9.98 x 5.08 mm<br>0.393 x 0.2" | 152.4 mm<br>6" |
|   | 25-1506 1PF DM |  | 4/50/1, 2 lbs          |                                   |                                |                |







## APPLICATOR STICKS

Made for general purpose functionality. Choose the size and shape that works best for your applications – specimen transfer, mixing rods, and laboratory work. Made with straight cut ends.

|   | Item No.  | Item Description | Packaging pack, weight | Dimensions length x diameter    |
|---|-----------|------------------|------------------------|---------------------------------|
|    | 805       | Wood.            | 12/1000, 23 lbs        | 304.8 x 2.108 mm<br>12 x 0.083" |
|   | 807       |                  | 20/1000, 18 lbs        |                                 |
|    | 25-807 2W |                  | 10/100/2, 6 lbs        |                                 |
|   | 809       | Wood.            | 12/500, 25 lbs         | 304.8 x 3.175 mm<br>12 x 0.125" |
|  | 811       | Wood.            | 2/1000, 12 lbs         | 457.2 x 3.175 mm<br>18 x 0.125" |

## CUTICLE/ORANGE STICKS

Single and double bevel ends offer numerous choices for all kinds of applications. Excellent for podiatric practices and long-term care. Ideal as an all-purpose stick for most clinical, diagnostic, and forensic settings. Products below shown at 40% of actual size.

|   | Item No. | Item Description                       | Packaging pack, weight | Dimensions length x diameter   |
|---|----------|--|------------------------|--------------------------------|
|  | 2900     | Wood, double bevel both ends.          | 5000 bulk, 16 lbs      | 177.8 x 3.962 mm<br>7 x 0.156" |
|   | 2910     |  | 10/100, 5 lbs          |                                |
|   | 2911     | Wood, double bevel both ends, wrapped. | 4/125/1, 6 lbs         |                                |
|  | 2920     | Wood, single bevel, double bevel ends. | 5000 bulk, 18 lbs      | 177.8 x 3.962 mm<br>7 x 0.156" |
|  | 2925     | Wood, single bevel, pointed one end.   | 10/100, 4 lbs          | 152.4 x 3.962 mm<br>6 x 0.156" |
|  | 2930     | Wood, single bevel, double bevel ends. | 10/100, 5 lbs          | 177.8 x 3.962 mm<br>7 x 0.156" |
|  | 2940     | Wood, double bevel both ends.          | 2500 bulk, 20 lbs      | 177.8 x 4.775 mm<br>7 x 0.188" |
|   | 2950     |  | 10/100, 6 lbs          |                                |
|  | 2970     | Wood, single bevel, double bevel ends. | 10/100, 6 lbs          | 177.8 x 4.775 mm<br>7 x 0.188" |

## TONGUE DEPRESSORS

While these are still called tongue depressors, wooden blades are used for a wide variety of applications, from patient oral exams to mixing or applying medicines, to crafts and industrial use.



Products below shown at 40% of actual size.

|  | Item No.       | Item Description                                 | Packaging pack, weight | Dimensions length x width         |
|--|----------------|--|------------------------|-----------------------------------|
|    | 700            | Infant.  | 10/1000, 37 lbs        | 114.3 x 9.53 mm<br>4.5 x 0.375"   |
|  | 25-700         |  | 10/100/1, 8 lbs        |                                   |
|    | 704            | Standard.  | 10/500, 35 lbs         | 152.4 x 17.475 mm<br>6 x 0.688"   |
|  | 704 BULK       |  | 3,700 bulk, 28 lbs     |                                   |
|    | 705            | Regular.   | 10/500, 38 lbs         | 152.4 x 19.05 mm<br>6 x 0.75"     |
|  | 25-705         |  | 10/100/1, 10 lbs       |                                   |
|    | 706            | Standard.  | 12/100, 9 lbs          | 152.4 x 17.475 mm<br>6 x 0.688"   |
|    | 709            | Standard, individually wrapped.                  | 10/250, 21 lbs         | 152.4 x 17.475 mm<br>6 x 0.688"   |
|  | 710            | Junior, individually wrapped, plain or flavored. | 10/250, 21 lbs         | 139.7 x 15.875 mm<br>5.5 x 0.625" |
|  | 710 CHERRY     |  | 20/50, 10 lbs          |                                   |
|  | 710 GRAPE      |  |                        |                                   |
|  | 710 STRAWBERRY |  |                        |                                   |
|  | 711            | Junior.  | 10/500, 29 lbs         | 139.7 x 15.875 mm<br>5.5 x 0.625" |
|  | 25-711         |  | 10/100/1, 10 lbs       |                                   |
|  | 712            |  | 12/100, 8 lbs          |                                   |

## MEDICAL SPOONS





Today Puritan Medical Products manufactures some of the most sophisticated single-use medical and diagnostic products on the market. But we still pay close attention to our simple wood products.



Products below shown at 70% of actual size.

|  | Item No. | Item Description        | Packaging pack, weight | Dimensions length x diameter   |
|--|----------|-------------------------|------------------------|--------------------------------|
|  | 625      | Wood spoon, double end. | 10/1000, 32 lbs        | 76.2 x 19.456 mm<br>3 x 0.766" |
|  | 645      | Wood spoon, double end. | 10/500, 31 lbs         | 127 x 25.4 mm<br>5 x 1"        |

## ECONOMY GRADE HOSPITAL & BRIGHTWOOD LINES

These lines offer lower-cost alternatives and include both US-made and imported materials. Spoons, blades, swabs, and applicators are ideal for mixing and dispensing medicines – handy for clinical kits and field work. Products below shown at 50% of actual size.

|   | Item No.                              | Item Description                                  | Packaging pack, weight | Dimensions length x width         |
|---|---------------------------------------|---|------------------------|-----------------------------------|
|  | 625 HOSPITAL<br>non-sterile   sterile | Wood spoon, double end.<br>Hospital® brand.       | 10/1000,<br>34 lbs     | 76.2 x 19.456 mm<br>3 x 0.766"    |
|  | 704 BRIGHTWOOD                        | Standard size tongue blade.<br>Brightwood® brand. | 10/500,<br>35 lbs      | 152.4 x 17.475 mm<br>6 x 0.688"   |
|  | 704 HOSPITAL                          | Standard size tongue blade.<br>Hospital® brand.   | 10/500,<br>35 lbs      | 152.4 x 17.475 mm<br>6 x 0.688"   |
|   | 25-705 HOSPITAL                       |   | 10/100/1,<br>11 lbs    |                                   |
|  | 711 HOSPITAL                          | Junior size tongue blade.<br>Hospital® brand.     | 10/500,<br>29 lbs      | 139.7 x 15.875 mm<br>5.5 x 0.625" |
|   | 25-711 HOSPITAL                       |   | 10/100/1,<br>9 lbs     |                                   |




|   | Item No.            | Item Description                                      | Packaging pack, weight | Handle length x diameter             | Tip length x diameter               | Full Length    |
|---|---------------------|---|------------------------|--------------------------------------|-------------------------------------|----------------|
|  | 803-WC HOSPITAL     | Standard cotton tip, wood handle.<br>Hospital® brand. | 10/10/100,<br>8 lbs    | 73.812 x 2.108 mm<br>2.906 x 0.083"  | 15.875 x 4.775 mm<br>0.625 x 0.188" | 76.2 mm<br>3"  |
|   | 25-803 2WC HOSPITAL |   | 10/100/1,<br>5 lbs     |                                      |                                     |                |
|  | 806-WC HOSPITAL     | Standard cotton tip, wood handle.<br>Hospital® brand. | 10/10/100,<br>14 lbs   | 147.624 x 2.108 mm<br>5.812 x 0.083" | 15.875 x 4.775 mm<br>0.625 x 0.188" | 152.4 mm<br>6" |
|   | 25-806 1WC HOSPITAL |   | 10/100/1,<br>6 lbs     |                                      |                                     |                |
|   | 25-806 2WC HOSPITAL |   | 10/100/2,<br>7 lbs     |                                      |                                     |                |

|   | Item No.     | Item Description                     | Packaging pack, weight | Dimensions length x diameter         |
|---|--------------|--------------------------------------|------------------------|--------------------------------------|
|  | 807 HOSPITAL | Wood applicator.<br>Hospital® brand. | 20/1000,<br>19 lbs     | 147.624 x 2.108 mm<br>5.812 x 0.083" |

## POPULES®

Super-convenient pre-filled popules with self-saturating, medical grade foam heads. Specify one of our standard solutions or let us know what you need.

Products below shown at 60% of actual size.

|   | Item No. | Item Description  | Packaging pack, weight | Handle length x diameter           | Tip length x diameter                | Full Length         |
|---|----------|---|------------------------|------------------------------------|--------------------------------------|---------------------|
|  | 4100     | Self-saturating foam swab, 3% hydrogen peroxide.                      | 4/3/50,<br>6 lbs       | 84.124 x 6.858 mm<br>3.312 x 0.27" | 28.575 x 15.875 mm<br>1.250 x 0.625" | 90.474 mm<br>3.562" |
|  | 4600 P   | Self-saturating foam swab, 10% povidone iodine prep solution.         | 4/150/1,<br>7 lbs      | 84.124 x 6.858 mm<br>3.312 x 0.27" | 34.925 x 17.449 mm<br>1.375 x 0.687" | 90.474 mm<br>3.562" |
|  | 4620     | Self-saturating foam swab, 70% isopropyl alcohol and deionized water. | 4/150/1,<br>7 lbs      | 84.124 x 6.858 mm<br>3.312 x 0.27" | 34.925 x 17.449 mm<br>1.375 x 0.687" | 90.474 mm<br>3.562" |





**Puritan**<sup>®</sup>  
Quality since 1919



**Puritan Medical**  
A Green Company

Puritan is passionate about being green.

We think green and give energy efficiency and resource stewardship high priority. Below are just a few things that make us green. Visit our website for more.

- Puritan runs on **100% renewable energy** (wood process waste) for all manufacturing energy needs.
- We have an **ongoing energy upgrade initiative**.
- Puritan packaging materials are made up of at least **65% post-consumer recycled content**.

**PURITAN MEDICAL PRODUCTS CO. LLC**  
P.O. BOX 149, 31 SCHOOL STREET  
GUILFORD, MAINE 04443-0149

Tel 800-321-2313 (US and Canada)  
207-876-3311

Fax 800-323-4153 (US and Canada)  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

#### THE PURITAN ADVANTAGE™

Every one of our products in the Puritan and PurSwab lines comes with the Puritan Advantage™ built right in.

We are a US manufacturer, with R&D, fabrication, sterilization, packaging, shipping, QC and oversight, sales, and customer service under one roof. Our products are consistently reliable and readily available.

We are a global company. Our medical devices carry the CE mark.

ISO 9001:2008 ISO 13485:2003 CE

#### GET SOCIAL WITH PURITAN!



[facebook.com/PuritanUSA](https://facebook.com/PuritanUSA)



[twitter.com/PuritanUSA](https://twitter.com/PuritanUSA)



[linkedin.com/company/PuritanUSA](https://linkedin.com/company/PuritanUSA)



[youtube.com/PuritanUSA](https://youtube.com/PuritanUSA)