

## MAPK7 Polyclonal Antibody

Catalog number: 10036-2-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

MAPK7; BMK1; ERK4; ERK5;

PRKM7



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 10036-2-AP (MAPK7 Antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse; other species not tested.
Positive WB detected in:	Mouse skeletal muscle tissue, HeLa cells, human heart tissue, mouse brain tissue
Calculated MAPK7 MW:	89kd
Observed MAPK7 MW:	59kd
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag0017
GenBank accession number:	BC007992
Gene ID (NCBI):	5598
Full name:	Mitogen-activated protein kinase 7

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.