

TIMM17A Polyclonal Antibody

Catalog number: 11189-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

TIMM17A; TIM17; TIM17A



Immunohistochemistry of paraffin-embedded human colon cancer using 11189-1-AP (TIMM17A Antibody) at Dilution 1:100 (under 10x lens)

Applications

Tested applications:	ELISA, IHC
Cited applications:	IHC, WB
Species specificity:	Human, mouse, rat; other species not tested.
Cited species:	Human
Calculated TIMM17A MW:	18kd
Observed TIMM17A MW:	
Positive IHC detected in:	Human colon cancer
Recommended dilution:	IHC: 1:50-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag1665
GenBank accession number:	BC015098
Gene ID (NCBI):	10440
Full name:	Translocase of inner mitochondrial membrane 17 homolog A (yeast)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.