

## PPP2R1B Polyclonal Antibody

Catalog number: 12621-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PPP2R1B; MGC26454; PR65B



A2780 cells were subjected to SDS PAGE followed by western blot with 12621-1-AP(PPP2R1B Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	A2780 cells, human kidney tissue, mouse brain tissue
Calculated PPP2R1B MW:	667aa,74kd
Observed PPP2R1B MW:	74kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag3307
GenBank accession number:	BC027596
Gene ID (NCBI):	5519
Full name:	Protein phosphatase 2 (formerly 2A), regulatory subunit A, beta isoform

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.