

CliniCool® Silver Series undercounter pharmacy/vaccine refrigerators are ideal for clinical labs, medical offices and pharmacies with a limited space and budget. These units offer superior temperature performance for safe storage of vaccines while meeting federal and state guidelines.



Standard Product Features:

- Cycle defrost
- Microprocessor temperature controller with battery back-up (batteries included)
- Adjustable operating temperature range: +36°F to +46°F (+2°C to +8°C)
- Freestanding thermometer with 2 year certificate of calibration
- Digital temperature display
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable
- Audible and visual high and low temperature alarms, with battery back up
- Remote alarm contacts
- 2 probes; 1 in air 1 in glycol bottle
- Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall
- Bottom footed shelf to allow air circulation
- Back wall evaporator cover
- Forced draft circulation
- Keyed door lock
- LED light
- Magnetic door gasket for positive seal
- Leveling legs
- Probe access port (3/8")
- Powder coat exterior
- High density urethane foam cabinet and door insulation
- HFC-free refrigerant (R600a), foam insulation and packaging
- Pharmacy refrigerator/freezer toolkit and temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage.
- Warranty: 2 years parts and labor



Hydrocarbon refrigerant used in all undercounter LHP models

Vaccine Compliant Features

All units tested and certified to maintain temperatures within refrigerated or frozen ranges required by the CDC

(1) Thermometer with 3 years NIST Certified of Calibration

Vaccine Storage Power Cord Warning Label

Vaccine Toolkit

Specifications and Ordering Information:

Model	LHP-2-UR-PH*	LHP-2-URG-PH*	LHP-3-UR-PH	LHP-3-URG-PH	LHP-3-URB-PH*	LHP-5-UR-PH
Weight lbs	63	63	68	68	69	93
Overall Height in.	21 1/4*	21 1/4*	28 7/8*	28 7/8*	30 3/4*	32 1/8*
Overall Width in.	17 1/4*	17 1/4*	17 3/4*	17 3/4*	17 3/4*	23 3/4*
Overall Depth in.	19 1/4*	19 1/4*	19 1/2*	19 1/2*	22 3/4*	24*
Interior cu. ft.	1	1	2.5	2.5	2.5	5.2
Number of Shelves	2, adj	2 adj	4 adj	4 adj	4 adj	3 adj
Door Type	solid	glass	solid	glass	solid	solid
Defrost Type	cycle	cycle	cycle	cycle	cycle	cycle
Compressor Size	1/8 HP	1/8 HP	1/8 HP	1/8 HP	1/8 HP	1/6HP
Plug (Receptacle)	5-15P	5-15P	5-15P	5-15P	5-15P	5-15P
Total Amp Draw	0.9	0.9	0.9	0.9	.9	1.3
Voltage	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz

*Left hinged model available

Model	LHP-5-URG-PH	LHP-5-URB-PH*	LHP-5-URBPH-ADA*	LHP-5-URBG-PH*	LHP-5-URBGPH-ADA*	LHP-5-URBSS-PH*
Weight lbs	93	97	97	97	97	100
Overall Height in.	32 1/8*	33 3/8*	31 15/16*	33 3/8*	31 15/16*	33 3/8*
Overall Width in.	23 3/4*	23 3/4*	23 3/4*	23 3/4*	23 3/4*	23 3/4*
Overall Depth in.	24*	24 1/2*	24 1/2*	25 1/2*	25 1/2*	24 1/2*
Interior cu. ft.	5.2	4.6	4.6	4.6	4.6	4.5
Number of Shelves	3 adj	2 adj, 1 fixed	2 adj, 1 fixed	2 adj, 1 fixed	2 adj, 1 fixed	2 adj, 1 fixed
Door Type	glass	solid	solid	glass	glass	Solid SS
Defrost Type	cycle	cycle	cycle	cycle	cycle	cycle
Compressor Size	1/6HP	1/8 HP	1/8 HP	1/8 HP	1/8 HP	1/8 HP
Plug (Receptacle)	5-15P	5-15P	5-15P	5-15P	5-15P	5-15P
Total Amp Draw	1.3	1.74	1.74	1.74	1.74	1.74
Voltage	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz

Exterior dimensions do not include 4" clearance required for installation.
Specifications subject to change without notice.



Receptacle Required



115V/1PH/ 60Hz
NEMA-5-15P



LABRepCo