

BUY 3 CASES OF KIMBLE® MEDIA BOTTLES AND GET 1 CASE FREE!*



GL 45 MEDIA BOTTLES

If you have never tried Kimble® KIMAX® GL 45 Media Bottles, now is the time. The bottles are autoclavable and ideal for general laboratory use including mixing, storing or transporting culture media, chemicals, & solvents.



***Free bottle capacity must be equivalent or lower than purchased models**

GL 45 Storage / Media Bottles With Blue Polypropylene Caps

KIMAX® bottles have a wide opening (approx. 30 mm ID) and GL 45 screw threads.

Each bottle is equipped with a clear drip-free, polypropylene pour ring and a liner-less, autoclavable (140°C max), polypropylene blue cap. KIMAX® bottles markings and graduations are made of permanent, white, enamel offering excellent resistance to heat, chemicals and solvents. These bottles are compatible with Kimble's broad range of caps and connection accessories for filtration, HPLC, etc...

Catalog Number	Capacity (mL)	Graduation Range (mL) x Interval	Approx. (Unassembled) OD - HT (mm)	Case Qty.
1743G50	100	20-80 x 10	56 - 100	10
1743G52	250	50-200 x 25	70 - 138	10
1743G54	500	100-400 x 50	86 - 176	10
1743G56	1000	100-900 x 50	101 - 225	10
1743G58	2000	400-1800 x 100	136 - 262	4
1743G57	5000	1000-4500 x 250	186 - 355	1
1743G59	10000	2000-9000 x 500	234 - 435	1

* Fill out the redemption form on the last page to receive your free product



We Believe You Are Important
— How Can We Help?

800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
thomassci.com • value@thomassci.com





Buy 3 Cases of Kimble® Media Bottles and get 1 Case Free!*

☐ **YES!** I have purchased 3 Cases of Kimble® Media Bottles.
Please send me 1 Free Case. Please select Thomas No. below:
**Free bottle capacity must be equivalent or lower than purchased models*

- | | |
|----------------------------------|----------------------------------|
| <input type="checkbox"/> 1743G50 | <input type="checkbox"/> 1743G58 |
| <input type="checkbox"/> 1743G52 | <input type="checkbox"/> 1743G57 |
| <input type="checkbox"/> 1743G54 | <input type="checkbox"/> 1743G59 |
| <input type="checkbox"/> 1743G56 | |

Offers subject to change without notice. Please allow 4-6 weeks for delivery of free products when applicable. Promotion expires December 31, 2016. All redemption forms must be received no later than January 14, 2017 to be valid.

Please complete the account information below,
and fax or email, along with Thomas Scientific Proof-of-Purchase, to:
Thomas Marketing: Fax 856-467-5167 Email Marketing@thomassci.com

Account Information

Name _____ Date _____
Organization _____
Account No. _____
Ship To Address _____
Room/Bldg./Dept. _____
City _____ State _____ Zip _____
Phone _____
Email _____



We Believe You Are Important
— How Can We Help?

800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
thomassci.com • value@thomassci.com

