

PAMCI Polyclonal Antibody

Catalog number: 12861-1-AP


Size: 40 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

RASSF9; P CIP1, PAMCI, PCIP1,
RASSF9



human liver tissue were subjected to SDS PAGE followed by western blot with 12861-1-AP(RASSF9 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated PAMCI MW:	435aa,50 kDa
Observed PAMCI MW:	50-55 kDa
Positive WB detected in	Human liver tissue, SKOV-3 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3976
GenBank accession number:	BC031589
Gene ID (NCBI):	9182
Full name:	Ras association (RalGDS/AF-6) domain family (N-terminal) member 9

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.