



Safely Fume Evidence with Consistent Results

Position your evidence inside the CA Chamber, add cyanoacrylate to the heating element and push the start button. Within minutes, high quality latent fingerprints are developed. Our exclusive, patented (# 6,660,054) design ensures potentially hazardous cyanoacrylate fumes generated during the development process never escape the chamber and must pass through the filtration system a minimum of twenty times before the purge process is complete. Precision humidity control up to 80%, a high temperature heating element and an internal circulation fan allow for complete control of the fuming process.



CA-3000

U.S. Patent No. 6,660,054

Extensive Safety Features:

- Patented, recirculatory filtration design
- Automatic locking door during development process
- Seamless door gasket for an air tight seal
- Special blended filter for cyanoacrylate fume capture
- Automatic monitoring of filter life — 500 cycles
- Automatic monitoring of pre-filter life — 50 cycles

Precision Microprocessor Controls:

- Fully automated development system
- Accurate control of relative humidity at 80%
- User definable fuming and filtering intervals
- LCD display of fuming options and system parameters
- Programmable internal circulation fan- promotes even cyanoacrylate fume distribution within chamber
- High temperature hot plate for rapid vaporization of cyanoacrylate

Design and Construction:

- Powder coated aluminum and steel super structure
- Safety glass viewing door
- Vapor proof internal fluorescent lightning
- Patented filtration design



CA-9000 High Capacity Fuming Chamber



CA-6000
(optional cart shown)

U.S. Patent No. 6,660,054

SPECIFICATIONS	CA-3000	CA-6000	CA-9000
Internal Dimensions	36"H x 30"W x 25"D	48"H x 30"W x 25"D	76"H x 70"W x 46"D
External Dimensions	53"H x 39"W x 32"D	65"H x 39"W x 32"D	94"H x 72"W x 48"D
Electrical Option	120V/ 60Hz or 230V/ 50Hz		
FILTERS:Units ship complete with filters; part numbers for reorder			
MY-2818 (CYN)	1 required	1 required	2 required
FE-2013 (Pre) pkg. of 6	1 required	1 required	2 required
CA-6000C Mobile steel cart	25"H x 39"W x 34"D		not applicable

* Extended warranty available, call for rates