


iClean Foam Swab, Foam Tip, 69mm Total Length

Thomas No.	Foam Tip		Total	
	Length	Diameter	Length	Breakpoint
20A00V257 Mfr. No. CY-FS742	12mm	3.5MM±0.2mm	69±2mm	Null



Product name	Foam swab
Swab Handle Material	Polypropylene
Swab Tip Material	polyurethane
Handle feature	no breakpoint but breakable
Connect Method	Hot Pressing
Function	Sample collection
Sterilization method	Irradiation

Thomas No. 20A00V257

iClean Foam Swab, Foam Tip, 69mm Total Length, 10/pack, 2000 swabs per case