

Low Lint Tissue Laboratory Wipers

- **LabForce® Low Lint Tissue Laboratory Wipers** are the perfect choice when your work requires delicate, light-duty wiping. Rely on LabForce® Wipers for cleaning surfaces, sensitive electronic parts, instruments in labs, laboratory lenses, and medical offices.
- The anti-static dispensing feature built-in to each package reduces electrostatic discharge which reduces dust particles remaining on the wiper, and allows for one-at-a-time dispensing.
- Perfect for the most sensitive jobs in the laboratory, the drafting room, the clinic, or any electronic applications
- Single-ply, pop-up box
- **LabForce® Low Lint Tissue Wipers** are available in various sizes and counts:
 - 280 ct Tissue Wipers measure 4.9" x 7.9" and is useful for cleaning laboratory instruments, polishing laboratory glassware and other delicate surfaces.
 - 140 ct Tissue Wipers measure 14.5" x 15.75" and is useful for larger cleaning areas.



Thomas No.	Description	Quantity
21A00A379	labForce Low Lint Delicate Task Wipes 4-7/8 x 7-7/8	280/Pk ; 16800/Cs
21A00A380	labForce Low Lint Delicate Task Wipes 14.5 x 16	140/Pk ; 2100/Cs