

CHRNA6 Polyclonal Antibody

Catalog number: 11388-1-AP


Size: 30 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CHRNA6; CHNRA6, CHRNA6



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 11388-1-AP(CHRNA6 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated CHRNA6 MW:	494aa,57 kDa
Observed CHRNA6 MW:	57 kDa
Positive WB detected in	Mouse thymus tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1921
GenBank accession number:	BC014456
Gene ID (NCBI):	8973
Full name:	Cholinergic receptor, nicotinic, alpha 6

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.