

Key Features:

- Manufactured using USP Class VI and USFDA 21 CFR Certified Material
- Leakproof when used with Foxx supplied caps
- Narrow Mouth for ease of use in automatic filling stations
- Wide Mouth for ease of sample collection and physical fill finish
- Suitable for Freeze Storage and General Storage Applications
- Liner Less Caps to reduce chances of secondary contamination
- The bottom notch ensures smooth labelling every time
- Low Liquid Retention of your precious reagent

Applications:

- Storage, Collection or Shipping of Chemicals, Buffers, Reagents or any liquid & powder samples
- Storage & Protection of Light-Sensitive Chemicals in Amber Bottle
- Wide Mouth Bottles for easy pouring and dispensing of viscous or powdered chemicals

The **EZLabpure™** Reagent Bottle series offers the perfect solution for safe, efficient, and reliable storage and transportation of reagents, buffer and chemical solutions for modern laboratories and testing kits. Designed for both liquids and powders, these bottles come in wide-mouth and narrow-mouth options to suit any handling need. All materials are USP Class VI-compliant and manufactured in a ISO Class 8 cleanroom, ensuring superior and consistent quality

Made from high-quality virgin PP, HDPE, and LDPE resins, they combine durability with excellent chemical resistance. PP bottles are autoclavable up to 121°C, ideal for sterilization before usage; HDPE bottles provide outstanding strength for long-term storage, while LDPE bottles are lightweight, flexible, and engineered for recurring shipments. The secure, high-grip cap ensures a completely leak-proof seal, while the collar of the bottle enables confident, stable handling in automatic filling stations even in extreme or gloved working conditions. Semi-buttress threads allow smooth, effortless opening and closing while maintaining a strong, and reliable closure fitment.

Available in sizes from 4 mL to 1000 mL, the bottles are versatile enough for any application, and their universal design makes them compatible with automated filling systems with high throughput.

Built to withstand extreme temperatures down to -100°C, the HDPE **EZLabpure™** Reagent Bottles deliver unmatched performance, reliability, and convenience, making them an essential choice for every lab and testing kits.

 **AUTOCLAVABLE**  **BPA-Free**




MADE  IN



EZLabpure™ Reagent Bottle Ordering Information:

Material	Part Number	Description	Volume	Qty/CS
HDPE	762-1111-FLS	EZLabpure™ HDPE Narrow Mouth Reagent Bottle, 4mL, Non-Sterile, 72/Case	4mL	72/Case
	762-1311-FLS	EZLabpure™ HDPE Narrow Mouth Reagent Bottle, 8mL, Non-Sterile, 72/Case	8mL	72/Case
	762-1511-FLS	EZLabpure™ HDPE Narrow Mouth Reagent Bottle, 15mL, Non-Sterile, 72/Case	15mL	72/Case
	762-1711-FLS	EZLabpure™ HDPE Narrow Mouth Reagent Bottle, 30mL, Non-Sterile, 72/Case	30mL	72/Case
	762-1811-FLS	EZLabpure™ HDPE Narrow Mouth Reagent Bottle, 60mL, Non-Sterile, 72/Case	60mL	72/Case
	762-1B11-FLS	EZLabpure™ HDPE Narrow Mouth Reagent Bottle, 125mL, Non-Sterile, 72/Case	125mL	72/Case
	762-1E11-FLS	EZLabpure™ HDPE Narrow Mouth Reagent Bottle, 250mL, Non-Sterile, 72/Case	250mL	72/Case
	762-1G21-FLS	EZLabpure™ HDPE Narrow Mouth Reagent Bottle, 500mL, Non-Sterile, 48/Case	500mL	48/Case
	762-1H31-FLS	EZLabpure™ HDPE Narrow Mouth Reagent Bottle, 1000mL, Non-Sterile, 24/Case	1000mL	24/Case
	762-2711-FLS	EZLabpure™ HDPE Wide Mouth Reagent Bottle, 30mL, Non-Sterile, 72/Case	30mL	72/Case
	762-2811-FLS	EZLabpure™ HDPE Wide Mouth Reagent Bottle, 60mL, Non-Sterile, 72/Case	60mL	72/Case
	762-2B11-FLS	EZLabpure™ HDPE Wide Mouth Reagent Bottle, 125mL, Non-Sterile, 72/Case	125mL	72/Case
	762-2E11-FLS	EZLabpure™ HDPE Wide Mouth Reagent Bottle, 250mL, Non-Sterile, 72/Case	250mL	72/Case
	762-2G21-FLS	EZLabpure™ HDPE Wide Mouth Reagent Bottle, 500mL, Non-Sterile, 48/Case	500mL	48/Case
	762-2H31-FLS	EZLabpure™ HDPE Wide Mouth Reagent Bottle, 1000mL, Non-Sterile, 24/Case	1000mL	24/Case
HDPE AMBER	76K-1111-FLS	EZLabpure™ HDPE Amber Narrow Mouth Reagent Bottle, 4mL, Non-Sterile, 72/Case	4mL	72/Case
	76K-1311-FLS	EZLabpure™ HDPE Amber Narrow Mouth Reagent Bottle, 8mL, Non-Sterile, 72/Case	8mL	72/Case
	76K-1511-FLS	EZLabpure™ HDPE Amber Narrow Mouth Reagent Bottle, 15mL, Non-Sterile, 72/Case	15mL	72/Case
	76K-1711-FLS	EZLabpure™ HDPE Amber Narrow Mouth Reagent Bottle, 30mL, Non-Sterile, 72/Case	30mL	72/Case
	76K-1811-FLS	EZLabpure™ HDPE Amber Narrow Mouth Reagent Bottle, 60mL, Non-Sterile, 72/Case	60mL	72/Case
	76K-1B11-FLS	EZLabpure™ HDPE Amber Narrow Mouth Reagent Bottle, 125mL, Non-Sterile, 72/Case	125mL	72/Case
	76K-1E11-FLS	EZLabpure™ HDPE Amber Narrow Mouth Reagent Bottle, 250mL, Non-Sterile, 72/Case	250mL	72/Case
	76K-1G21-FLS	EZLabpure™ HDPE Amber Narrow Mouth Reagent Bottle, 500mL, Non-Sterile, 48/Case	500mL	48/Case
	76K-1H31-FLS	EZLabpure™ HDPE Amber Narrow Mouth Reagent Bottle, 1000mL, Non-Sterile, 24/Case	1000mL	24/Case
	76K-2711-FLS	EZLabpure™ HDPE Amber Wide Mouth Reagent Bottle, 30mL, Non-Sterile, 72/Case	30mL	72/Case
	76K-2811-FLS	EZLabpure™ HDPE Amber Wide Mouth Reagent Bottle, 60mL, Non-Sterile, 72/Case	60mL	72/Case
	76K-2B11-FLS	EZLabpure™ HDPE Amber Wide Mouth Reagent Bottle, 125mL, Non-Sterile, 72/Case	125mL	72/Case
	76K-2E11-FLS	EZLabpure™ HDPE Amber Wide Mouth Reagent Bottle, 250mL, Non-Sterile, 72/Case	250mL	72/Case
	76K-2G21-FLS	EZLabpure™ HDPE Amber Wide Mouth Reagent Bottle, 500mL, Non-Sterile, 48/Case	500mL	48/Case
	76K-2H31-FLS	EZLabpure™ HDPE Amber Wide Mouth Reagent Bottle, 1000mL, Non-Sterile, 24/Case	1000mL	24/Case



Material	Part Number	Description	Volume	Qty/CS
PP	765-1111-FLS	EZLabpure™ PP Narrow Mouth Reagent Bottle, 4mL, Non-Sterile, 72/Case	4mL	72/Case
	765-1311-FLS	EZLabpure™ PP Narrow Mouth Reagent Bottle, 8mL, Non-Sterile, 72/Case	8mL	72/Case
	765-1511-FLS	EZLabpure™ PP Narrow Mouth Reagent Bottle, 15mL, Non-Sterile, 72/Case	15mL	72/Case
	765-1711-FLS	EZLabpure™ PP Narrow Mouth Reagent Bottle, 30mL, Non-Sterile, 72/Case	30mL	72/Case
	765-1811-FLS	EZLabpure™ PP Narrow Mouth Reagent Bottle, 60mL, Non-Sterile, 72/Case	60mL	72/Case
	765-1B11-FLS	EZLabpure™ PP Narrow Mouth Reagent Bottle, 125mL, Non-Sterile, 72/Case	125mL	72/Case
	765-1E11-FLS	EZLabpure™ PP Narrow Mouth Reagent Bottle, 250mL, Non-Sterile, 72/Case	250mL	72/Case
	765-1G21-FLS	EZLabpure™ PP Narrow Mouth Reagent Bottle, 500mL, Non-Sterile, 48/Case	500mL	48/Case
	765-1H31-FLS	EZLabpure™ PP Narrow Mouth Reagent Bottle, 1000mL, Non-Sterile, 24/Case	1000mL	24/Case
	765-2711-FLS	EZLabpure™ PP Wide Mouth Reagent Bottle, 30mL, Non-Sterile, 72/Case	30mL	72/Case
	765-2811-FLS	EZLabpure™ PP Wide Mouth Reagent Bottle, 60mL, Non-Sterile, 72/Case	60mL	72/Case
	765-2B11-FLS	EZLabpure™ PP Wide Mouth Reagent Bottle, 125mL, Non-Sterile, 72/Case	125mL	72/Case
	765-2E11-FLS	EZLabpure™ PP Wide Mouth Reagent Bottle, 250mL, Non-Sterile, 72/Case	250mL	72/Case
	765-2G21-FLS	EZLabpure™ PP Wide Mouth Reagent Bottle, 500mL, Non-Sterile, 48/Case	500mL	48/Case
	765-2H31-FLS	EZLabpure™ PP Wide Mouth Reagent Bottle, 1000mL, Non-Sterile, 24/Case	1000mL	24/Case
LDPE	764-1311-FLS	EZLabpure™ LDPE Narrow Mouth Reagent Bottle, 8mL, Non-Sterile, 72/Case	8mL	72/Case
	764-1511-FLS	EZLabpure™ LDPE Narrow Mouth Reagent Bottle, 15mL, Non-Sterile, 72/Case	15mL	72/Case
	764-1711-FLS	EZLabpure™ LDPE Narrow Mouth Reagent Bottle, 30mL, Non-Sterile, 72/Case	30mL	72/Case
	764-1811-FLS	EZLabpure™ LDPE Narrow Mouth Reagent Bottle, 60mL, Non-Sterile, 72/Case	60mL	72/Case
	764-1B11-FLS	EZLabpure™ LDPE Narrow Mouth Reagent Bottle, 125mL, Non-Sterile, 72/Case	125mL	72/Case
	764-1E11-FLS	EZLabpure™ LDPE Narrow Mouth Reagent Bottle, 250mL, Non-Sterile, 72/Case	250mL	72/Case
	764-1G21-FLS	EZLabpure™ LDPE Narrow Mouth Reagent Bottle, 500mL, Non-Sterile, 48/Case	500mL	48/Case
	764-1H31-FLS	EZLabpure™ LDPE Narrow Mouth Reagent Bottle, 1000mL, Non-Sterile, 24/Case	1000mL	24/Case
	764-2711-FLS	EZLabpure™ LDPE Wide Mouth Reagent Bottle, 30mL, Non-Sterile, 72/Case	30mL	72/Case
	764-2811-FLS	EZLabpure™ LDPE Wide Mouth Reagent Bottle, 60mL, Non-Sterile, 72/Case	60mL	72/Case
	764-2B11-FLS	EZLabpure™ LDPE Wide Mouth Reagent Bottle, 125mL, Non-Sterile, 72/Case	125mL	72/Case
	764-2E11-FLS	EZLabpure™ LDPE Wide Mouth Reagent Bottle, 250mL, Non-Sterile, 72/Case	250mL	72/Case
	764-2G21-FLS	EZLabpure™ LDPE Wide Mouth Reagent Bottle, 500mL, Non-Sterile, 48/Case	500mL	48/Case
	764-2H31-FLS	EZLabpure™ LDPE Wide Mouth Reagent Bottle, 1000mL, Non-Sterile, 24/Case	1000mL	24/Case