

PNKP Polyclonal Antibody

Catalog number: 13558-1-AP


Size: 35 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PNKP; DNA 5 kinase/3
phosphatase, PNK, PNKP



Jurkat cells were subjected to SDS
PAGE followed by western blot with
13558-1-AP(PNKP Antibody) at
dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated PNKP MW:	521aa,57 kDa
Observed PNKP MW:	57 kDa
Positive WB detected in	Jurkat cells, HeLa cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC= Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4421
GenBank accession number:	BC033822
Gene ID (NCBI):	11284
Full name:	Polynucleotide kinase 3'-phosphatase

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.