

IN STOCK! Ready to Ship!



Cardboard Freezer Boxes with Dividers

labForce cryo/freezer boxes and dividers are made of sturdy white cardboard, with a moisture-repellent coating to protect from the effects of liquid nitrogen or condensation from refrigerated equipment. Boxes are sized to fit standard stainless-steel and polycarbonate storage systems. Small holes on the bottom allow for safe drainage of liquid nitrogen or water. For use in mechanical freezers and liquid nitrogen tanks at temperature to -196°C.

Product Features:

- With a moisture-repellent coating to protect from condensation
- A grid on the inside of the box (81- and 100-place boxes) for easy lining up and location tracking

Thomas No.	Description	Quantity
1141Q92	labForce 2" Cardboard Freezer Box with 81-Place Cell Divider	10/Pack
		50/Case
1141Q93	labForce 2" Cardboard Freezer Box with 100-Place Cell Divider	10/Pack
		50/Case