

Product Name:	ABS SSC 20
Description:	<p>ABS SSC Class Includes 6 Canisters Total 2ml Vial Capacity: 210 Liquid Nitrogen Capacity: 20.5 Liters ½ CC Straw (1 Level Bulk): 1122 ½ CC Straw (10 Cane): 720 Static Evaporation Rate (Liters/day): 0.09 Normal Working Time: 135 Days Neck Diameter (mm): 55.4 Exterior Height (mm): 660 Exterior Diameter (mm): 368.3 Canister Diameter (mm): 41.9 Durable aluminum tank with powder coated exterior Lid lock provision fiberglass restricted opening neck design Compatible design for both liquid and vapor phase storage LN2 measuring rod for level validation Color-coded canisters for inventory identification One-year parts and labor warranty Three-year vacuum warranty</p>

