



Item No.	Qty	Description	Inventory Code	Material
D1	1	50mL SU Storage Bags	ACAX0012	ULDPE
B1	2	1/8" 1/4" (3.2*6.4mm) L-Flex Tubing, Length 300mm	DBBX0001	TPE
C1	1	1/8" (3.2mm) Female Luer Thread Barb	EAAX0140	PP
C2	1	Male luer integral lock ring plug	EAAX0137	PP
C3	1	1/8" (3.2mm) Male luer integral barb	EAAX0132	PP
C4	1	Female luer thread style cap	EAAX0221	PP
F1	4	2.5mm cable tie	ADAA0005	Nylon
F2	2	Robert clip (small)	EDBA0026	PA

<b>Customer Approval</b>		<b>LePure Biotech Co., Ltd.</b>				
Company:		Drawing No. <b>HP0050-1AK</b>				
Signature:		Title: <b>SU Storage Bags</b>				
Date:		<b>Designed By:</b>	<b>Date:</b>	<b>Checked By:</b>	<b>Date:</b>	<b>Edition</b>
		Jiajia Guo	2023.03.12	Alan	2023.03.12	Rev01
<small>Design Approval the manufacturing of the system(s) based on this design. By signing, you also assume full</small>		<small>We request that you do not disclose this information to any third party involved in this drawing, in order to protect our customer's confidentiality and intellectual property</small>				