



## Aseptic Packing Line Medium

M1350

Aseptic Packing Line Medium is a fluid medium used in validating aseptic packing lines.

### Composition\*\*

Ingredients	Gms / Litre
Peptic digest of animal tissue	5.000
Yeast extract	2.500
Sucrose	5.000
Sodium chloride	5.000
Phenol red	0.005
Final pH ( at 25°C)	7.3±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 17.50 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense in tubes containing inverted Durhams tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

### Principle And Interpretation

Aseptic Packing Line Medium is used in validating aseptic packing and filling lines. It is a basic medium for demonstrating growth by either acid or gas production.

Peptic digest of animal tissue and yeast extract provide basic nutrients. Sucrose is the carbohydrate source. Sodium chloride maintains osmotic balance. Phenol red is the pH indicator. Sucrose fermentation leads to the formation of acidic conditions which is detected by phenol red, observed as a colour change from red to yellow.

- Dispense the reconstituted medium into packing line upstream of the sterilization process.
- Incubate the filled packs at 30±2°C for upto 7 days.

Gas production is observed as trapped air bubbles in the Durhams tubes and acid production is visualized as a colour change of the medium to yellow. Growth is indicated as turbidity in the medium.

### Quality Control

#### Appearance

Light yellow to light pink homogeneous free flowing powder

#### Colour and Clarity of prepared medium

Orange red coloured clear solution in tubes

#### Reaction

Reaction of 1.75% w/v aqueous solution at 25°C. pH : 7.3±0.2

#### pH

7.10-7.50

#### Cultural Response

M1350: Cultural characteristics observed after an incubation at 30±2°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth	Acid	Gas
<b>Cultural Response</b> <i>Bacillus cereus</i> ATCC 14579 50-100		luxuriant	positive reaction, yellow colour	negative reaction

<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	positive reaction, yellow colour	positive reaction
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	positive reaction, yellow colour	positive reaction
<i>Staphylococcus aureus</i> ATCC 25923	50-100	luxuriant	positive reaction, yellow colour	negative reaction

### Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on label.

Revision : 1 / 2011



#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.