

AK5 Polyclonal Antibody

Catalog number: 12510-2-AP


Size: 39 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

AK5; AK6; MGC33326



human brain tissue were subjected to SDS PAGE followed by western blot with 12510-2-AP(AK5 Antibody) at dilution of 1:2000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human brain tissue
Calculated AK5 MW:	63kd
Observed AK5 MW:	63kd
Recommended dilution:	WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3203
GenBank accession number:	BC012467
Gene ID (NCBI):	26289
Full name:	Adenylate kinase 5

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.