

GNL2 Polyclonal Antibody

Catalog number: 10190-2-AP

Size: 60 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

GNL2; FLJ40906;

HUMAQUANTIG; NGP1; Ngp-1



Immunohistochemistry of paraffin-embedded human breast cancer using 10190-2-AP (GNL2 Antibody) at Dilution 1:50 (under 10x lens)

Applications

Tested applications:	ELISA, IHC
Species specificity:	Human; other species not tested.
Calculated GNL2 MW:	731aa,84kd
Observed GNL2 MW:	
Positive IHC detected in:	Human breast cancer
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag0243
GenBank accession number:	BC000107
Gene ID (NCBI):	29889
Full name:	Guanine nucleotide binding protein-like 2 (nucleolar)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.