

DNAJA2 Polyclonal Antibody

Catalog number: 12236-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

DNAJA2; DNAJA2, DJ3, DJA2,

DNAJ, PRO3015



A375 cells were subjected to SDS PAGE followed by western blot with 12236-1-AP (DNAJA2 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Cited applications:	WB
Species specificity:	Human, mouse, rat; other species not tested.
Cited species:	Human
Calculated DNAJA2 MW:	417aa, 46 kDa
Observed DNAJA2 MW:	46kd
Positive WB detected in	A375 cells, human brain tissue, PC-3 cells
Recommended dilution:	WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2877
GenBank accession number:	BC015809
Gene ID (NCBI):	10294
Full name:	DnaJ (Hsp40) homolog, subfamily A, member 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.