

Cadherin-18 Polyclonal Antibody

Catalog number: 13091-1-AP

Size: 58 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

CDH18; cadherin 18, type 2,

Cadherin-14, Cadherin-18,

CDH14, CDH14L, CDH18,

CDH24, EY CADHERIN



COLO 320 cells were subjected to SDS PAGE followed by western blot with 13091-1-AP (CDH18 antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated Cadherin-18 MW:	790aa, 88kd
Observed Cadherin-18 MW:	105kd
Positive WB detected in	COLO 320 cells, A549 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3675
GenBank accession number:	BC031051
Gene ID (NCBI):	1016
Full name:	Cadherin 18, type 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.