

SLAMF7 Polyclonal Antibody

Catalog number: 12905-1-AP


Size: 38 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

SLAMF7; 19A, CD2 subset 1,
CD319, CRACC, CS1, Membrane
protein FOAP 12, Novel Ly9,
Protein 19A, SLAM family
member 7, SLAMF7



Jurkat cells were subjected to SDS
PAGE followed by western blot with
12905-1-AP(SLAMF7 antibody) at
dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated SLAMF7 MW:	335aa,37 kDa
Observed SLAMF7 MW:	37 kDa
Positive WB detected in	Jurkat cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3804
GenBank accession number:	BC027867
Gene ID (NCBI):	57823
Full name:	SLAM family member 7

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.