

TrkC Polyclonal Antibody

Catalog number: 11999-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

NTRK3; GP145 TrkC,
gp145(trkC), NT 3 growth factor
receptor, NTRK3, TrkC, TrkC
tyrosine kinase



human brain tissue were subjected to SDS PAGE followed by western blot with 11999-1-AP(NTRK3 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated TrkC MW:	839aa,94kd
Observed TrkC MW:	68kd and 94kd and 145kd
Positive WB detected in	Human brain tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2635
GenBank accession number:	BC013693
Gene ID (NCBI):	4916
Full name:	Neurotrophic tyrosine kinase, receptor, type 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.