


GEB Robotic Tips, compatible with Rainin® LTS® pipettes, can be used in proteomics, drug development, genomics and other fields.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (4.5*3.5)mm	Packaging Spec.
AT200R-B-N	200ul	Natural				1000tips/bag, 20 bags/case, 20000 tips/case.
AT200R-B-N-LB	200ul	Natural	Low Retention			
Case dimension L*W*H: 595*325*310mm			Gross Weight: 3.2 kg			
AF200R-B-N	200ul	Natural			Filtered	1000 tips/bag, 10 bags/case, 10000 tips/case.
AF200R-B-N-LB	200ul	Natural	Low Retention		Filtered	
Case dimension L*W*H:480* 345*280mm			Gross Weight: 2.0 kg			
AT200R-E-N	200ul	Natural				96 tips/rack,10 stacks /bag, 10bags/case,9600tips/case.
AT200R-E-N-LB	200ul	Natural	Low Retention			
Case dimension L*W*H:455*270 *265mm			Gross Weight: 6.3 kg			
AT200R-9-N	200ul	Natural				96 tips/rack, 10 racks /pack, 5packs/ case.50 racks/case, 4800tips/case.
AT200R-9-N-LB	200ul	Natural	Low Retention			
AT200R-9-NS	200ul	Natural		Sterilized		
AT200R-9-NS-LB	200ul	Natural	Low Retention	Sterilized		
AF200R-9-NS	200ul	Natural		Sterilized	Filtered	
AF200R-9-N-LB	200ul	Natural	Low Retention		Filtered	
AF200R-9-NS-LB	200ul	Natural	Low Retention	Sterilized	Filtered	
AT200R-9B-N	200ul	Natural				
AT200R-9B-NS	200ul	Natural		Sterilized		
AT200R-9B-N-LB	200ul	Natural	Low Retention			
AT200R-9B-NS-LB	200ul	Natural	Low Retention	Sterilized		
AF200R-9B-N	200ul	Natural			Filtered	
AF200R-9B-N-LB	200ul	Natural	Low Retention		Filtered	
AF200R-9B-NS	200ul	Natural		Sterilized	Filtered	
AF200R-9B-NS-LB	200ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H:480* 345*280mm			Gross Weight: 7.8kg			
AF200R-R-NS	200ul	Natural		Sterilized	Filtered	96 tips/rack, 10 racks /pack, 5 packs / case, 50 racks/case 4800tips/case.
AF200R-R-NS-LB	200ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H:500*265 *450mm			Gross Weight: 9.1kg			

- ◇ 200μL Natural Tips
- ◇ 96 Tip/Hanging Tip Rack format/Blut
- ◇ Polypropylene
- ◇ Certified DNase/RNase free
- ◇ Nonpyrogenic
- ◇ Presterile/Nonsterile
- ◇ Filtered (Filtered With aerosol-resistant filters) /NonFiltered
- ◇ Endotoxin-Safe meets USP/FDA standards
- ◇ Comply With USP 3127 (Heavy Metal Contamination) standards
- ◇ Compatible with Rainin® LTS® pipettes.
- ◇

Drawing: Reference tolerance table for unmarked tolerance



Select	± A	± B	± C
0-6	0.02	0.05	0.10
>6-30	0.05	0.10	0.15
>30-60	0.10	0.15	0.20
>60-120	0.15	0.20	0.25
>120-180	0.20	0.25	0.30
>180-300	0.25	0.30	0.35

Prepared by: Qun Zhang 2021/8/12

Approve: WeiMing Yu 2021/8/12