

OmniGenX Microplate Seal Pads



***The Fastest, Most
Inexpensive Sealing
Solution on the Market!***

Features:

- Inexpensive alternative to standard plate seals
- Easy to use for Fast and Effective plate sealing
- Ideal for sealing ELISA and Filtration plates
- Prevent Cross-Contamination and Evaporation between Aspiration and Dispense Steps

The OmniGenX Tape Pad is an inexpensive alternative to standard adhesive plate seals and lids. It's very easy to use. Simply remove a sheet from the pad and apply it to your plate or block. **There is NO backing to remove!!**

These seals are a great low cost solution for sealing microplates, ELISA and filtration plates, and for preventing cross-contamination and evaporation between aspiration and dispense steps. They are NOT designed for PCR or cold storage applications.

Thomas No.	Description	Qty.
1149J29	OmniGenX Microplate Seal Pad; 25 Seals/Pad, 4 Pads/Pack	100/PK 5000/Cs

Visit our website to watch a [demonstration video](#) and to place your order!