

# RHINOSWAB JUNIOR SPECIFICATIONS

RHINOSWAB  
JUNIOR  
by RHINOMED

**Rhinoswab Junior** delivers comfort, reliability and performance, enhanced by novel features that reduce the fear and distress associated with the use of respiratory swabs for children aged 4 years and over.

Rhinoswab Junior can be used for rapid antigen and PCR tests, working seamlessly with existing vials and transport media.

**Studies have confirmed that the Rhinoswab Junior<sup>1</sup>**

1. Has comparable sensitivity to regular swabs across 12 different respiratory targets.
2. Is preferred by children and their parents over regular swabs.

## RHINOSWAB JUNIOR REGULATORY INFORMATION:

**FDA Registration Number:** 3011194487

**ARTG Identifier:** 349491

**MHRA Application Reference:**

**Health Canada MDEL:** 6789

**CE Mark Certificate Number:** 542839

2022042501259656

### APPLICATIONS

- For collection of pathogens and biological samples from the nose for diagnostic testing of respiratory viruses including COVID-19
- Designed for easy self collection for rapid antigen and PCR tests

### DIMENSIONS

- 44.3mm total length
- 34.0mm total width
- 4.5mm nasal tip width
- 6.4mm nasal tip height
- 16.0mm total width of double loop

### QUALITIES

- Larger contact area to maximize sample collection
- Superior sample elution
- Medical grade quality
- Painless and simple self-administration

**Loops** ergonomically designed to sit comfortably in each nostril for pain free self sampling

**Loops** break off into standard collection tubes

**Double loops** anatomically designed to maximise contact area between the swab and the nose to maximise yield

**Novelty handle** allows for easy and standardised self-swabbing and reduces testing fears



Thomas No. 23A00A011



Thomas No. 23A00A008



Thomas No. 23A00A009



Thomas No. 23A00A010

**CONTACT US TO DISCUSS HOW RHINOSWAB MAY HELP DELIVER A  
BETTER SAMPLING SOLUTION FOR YOUR NEEDS**

<sup>1</sup>Tosif S, Lee L, Nguyen J, Overmars I, Selman C, Grobler A, McMinn A, Waller G, McNab S, Jarvis T, Steer A, Babi F, Daley A, Crawford N, A novel anterior nasal swab to detect respiratory viruses. In: Communicable Diseases & Immunisation Conference 2022; 20-22nd June, Sydney.