

ACVR1B Polyclonal Antibody

Catalog number: 10086-1-AP


Size: 50 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

ACVR1B; activin A receptor,
type IB, Activin receptor like
kinase 4, Activin receptor type
1B, Activin receptor type IB,
ACTR IB, ACTR1B, ACVR1B,
ACVRLK4, ALK 4, ALK4, SKR2



HEK-293 cells were subjected to SDS PAGE followed by western blot with 10086-1-AP (ACVR1B Antibody) at dilution of 1:300

Applications

| | |
|-------------------------|--|
| Tested applications: | ELISA, WB |
| Species specificity: | Human, Mouse, Rat; other species not tested. |
| Calculated ACVR1B MW: | 57 kDa |
| Observed ACVR1B MW: | 49 kDa |
| Positive WB detected in | HEK-293 cells, DMSO treated HepG2 cells |
| Recommended dilution: | WB: 1:200-1:1000 |

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

| | |
|---------------------------|-----------------------------|
| Immunogen: | Ag0129 |
| GenBank accession number: | BC000254 |
| Gene ID (NCBI): | 91 |
| Full name: | Activin A receptor, type IB |

Product information

| | |
|----------------------|--|
| Purification method: | Antigen affinity purification |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C. |