

SQUARE MEDIA BOTTLES

BUY 2 GET

1 FREE!



Square bottle shape conserves valuable space in refrigerated storage areas

- PETG materials reduce gas permeability (compared to polystyrene) to improve solution stability
- Lightweight plastic material is a safer, shatterproof alternative to glass media bottles
- Ideal for storing media, reagents and buffers

Part No.	Product Description	Qty/Case
Square Media Bottles		
24A00A639	30mL Media Bottle, Square	24
24A00A640	60mL Media Bottle, Square	24
24A00A641	125mL Media Bottle, Square	24
24A00A642	250mL Media Bottle, Square	24
24A00A643	500mL Media Bottle, Square	24
24A00A644	1000mL Media Bottle, Square	24
230756	2L Media Bottle, Square	12
Square Media Bottles — Tray Packed		
24A00A645	30mL Media Bottle, Square	48
24A00A646	60mL Media Bottle, Square	48
24A00A647	125mL Media Bottle, Square	24
24A00A648	250mL Media Bottle, Square	24
24A00A649	500mL Media Bottle, Square	24
24A00A650	1000mL Media Bottle, Square	24
230766	2L Media Bottle, Square	12
Square Media Bottles — Tray Packed, Tamper Evident		
25A00A210	30mL Media Bottle, Square	48
25A00A211	60mL Media Bottle, Square	48
25A00A212	125mL Media Bottle, Square	24
25A00A213	250mL Media Bottle, Square	24
25A00A214	500mL Media Bottle, Square	24
25A00A215	1000mL Media Bottle, Square	24
230766TE	2L Media Bottle, Square	12

Buy 2 of same product to get 1 Free of your choice

How To Redeem:

OPTION 1:



Submit online at:
celltreat.us/submitpromo

PROMO CODE:

B2G1-SQUARE-TSI26-Q3

Offer valid in the U.S. only and void where prohibited by government, law or company policy. Offer valid for end-user customers only for orders placed between July 1, 2026 and September 30, 2026. Limit 4 free cases per lab. Allow 4 - 6 weeks for delivery of product. Offer does not apply to previous purchases and cannot be combined with any other CELLTREAT/VistaLab discounts or promotional offers. Submission form must be received by October 15, 2026.



833.544.SHIP (7447)
CustomerService@ThomasSci.com
www.ThomasSci.com

OPTION 2:

Email form & proof of purchase to
promotions@celltreat.com

Name: _____

Company: _____

Shipping Address: _____

Lab/Room #: _____

City: _____ State: _____ Zip: _____

Email: _____

FREE CASE SELECTION - Eligible part numbers above, limit 4 FREE cases per lab

Case 1		Case 3	
Case 2		Case 4	