

PRRG2 Polyclonal Antibody

Catalog number: 11618-1-AP

Size: 56 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

PRRG2; PRGP2, Proline rich

Gla protein 2, PRRG2, TMG2



A549 cells were subjected to SDS PAGE followed by western blot with 11618-1-AP (PRRG2 antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated PRRG2 MW:	179aa, 20 kDa
Observed PRRG2 MW:	16 kDa
Positive WB detected in:	A549 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2197
GenBank accession number:	BC026032
Gene ID (NCBI):	5639
Full name:	Proline rich Gla (G-carboxyglutamic acid) 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.