

MOBP Polyclonal Antibody

Catalog number: 12472-1-AP


Size: 22 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

MOBP; MGC87379



mouse brain tissue were subjected to SDS PAGE followed by western blot with 12472-1-AP(MOBP Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse brain tissue
Calculated MOBP MW:	183aa,21kd
Observed MOBP MW:	10kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3156
GenBank accession number:	BC022471
Gene ID (NCBI):	4336
Full name:	Myelin-associated oligodendrocyte basic protein

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.