

## STAT4 Polyclonal Antibody

Catalog number: 13028-1-AP

Size: 24 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

STAT4; SLEB11



Jurkat cells were subjected to SDS PAGE followed by western blot with 13028-1-AP (STAT4 Antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Jurkat cells, HeLa cells, PC-3 cells
Calculated STAT4 MW:	748aa, 86kd
Observed STAT4 MW:	86kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag3660
GenBank accession number:	BC031212
Gene ID (NCBI):	6775
Full name:	Signal transducer and activator of transcription 4

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.