

ZC3HC1 Polyclonal Antibody

Catalog number: 13148-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

ZC3HC1; 1110054L24Rik;

HSPC216; NIPA



mouse lung tissue were subjected to SDS PAGE followed by western blot with 13148-1-AP(ZC3HC1 Antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse lung tissue
Calculated ZC3HC1 MW:	502aa,55kd
Observed ZC3HC1 MW:	55kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3744
GenBank accession number:	BC028917
Gene ID (NCBI):	51530
Full name:	Zinc finger, C3HC-type containing 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.